

GUANTES DE LATEX CON POLVO AQL 1.5

(Ligeramente empolvado)



TALLA	CODIGO	EAN-13
P	066093	8 435125 401890
M	066094	8 435125 401906
G	066095	8 435125 401913

CARACTERISTICAS

- LATEX HIPOALERGENICO
- ANATOMICO
- ANTIDESLIZANTE
- RESISTENTE Y CON SENSIBILIDAD AL TACTO
- CON REBORDE
- FORMA PLANA, CON SUPERFICIE EXTERNA LISA
- FABRICADOS BAJO NORMA UNE 455/1/2/3 – EN 420 y EN 374-2
- MARCA: STAR/PREMIUM
- FABRICANTE: DRGN LATEX GROUP SDN. BHD.
- PROCEDENCIA: MALAYSIA
- MARCADO “CE” : CLASE I , DECLARACION DE CONFORMIDAD

(1) Dimensiones

UNE 455 /1/2/3

Talla super-pequeña

1.- Longitud (mm.)	240 minimo / 260 mm.	240 min +/- 10 mm.
2.- Anchura (mm.)	78 + / - 3	78 + / -10
3.- Dedo (mm.)	0.25 + / -0.03	0.08 min
4.- Palma (mm.)	0.25 + / -0.03	0.08 min

Talla pequeña

1.- Longitud (mm.)	240 minimo / 260 mm.	240 min +/- 10 mm.
2.- Anchura (mm.)	82 + / - 3	80 + / -10
3.- Dedo (mm.)	0.25 + / -0.03	0.08 min
4.- Palma (mm.)	0.25 + / -0.03	0.08 min

Talla mediana

1.- Longitud (mm.)	240 minimo / 260 mm.	240 min+ /- 10 mm.
2.- Anchura (mm.)	97 + / - 3	95 + / -10
3.- Dedo (mm.)	0.25 + / -0.03	0.08 min
4.- Palma (mm.)	0.25 + / -0.03	0.08 min

(1) Dimensiones**UNE 455 /1/2/3****Talla grande**

1.- Longitud (mm.)	240 minimo / 260 mm.	240 min +/- 10 mm.
2.- Anchura (mm.)	112 + / - 3	110 + / -10 mm.
3.- Dedo (mm.)	0.25 + / -0.03	0.08 min
4.- Palma (mm.)	0.25 + / -0.03	0.08 min

(2) Propiedades de Tensión**UNE 455 /1/2/3**

Resistencia a la tensión (Mpa.)	25 + / - 3.0	23 min
Elongación (en %)	820 + / - 30	500 a 700 máximo
Lubricante	Polvos de arroz o Maiz	Polvos de arroz o maiz. Nunca polvos de talco.
Modulus 500 % (Mpa.)	Máximo 5.5	
Fuerza a la rotura (N)	Mínimo 10,5	10,5 Newtons

(3) Diseño**UNE 455 /1/2/3**

Ambidextros.	Ambidextros.
--------------	--------------

(4) Nivel de proteínas

Confirmado, menor 150 µg./gr. Normalmente, entre 70 y 110 µg./gr.
--

(5) Ligeramente empolvado con almidón de maiz**PRESENTACION**

	UNE 455 /1/2/3
Caja dispensadora 100 Uds. 10 x 100 Uds. / Cartón	Caja dispensadora 100 Uds.

ALMACENAJE: A TEMPERATURA AMBIENTE

DRGN LATEX GROUP SDN. BHD.

N° 46-1, Jalan Nautika BU20/B, Pusat Komersial TSB.

Sungai Buloh, 47000 Selangor, Malaysia

Phone 03 – 339 31563

LATEX EXAMINATION GLOVES, POWDERED, NON-STERILE

DECLARATION OF CONFORMITY

LATEX EXAMINATION GLOVES POWDERED

We, like manufacturer of **LATEX EXAMINATION GLOVES POWDERED** and according to European Council Directive 93/42/EEC dated on 14 June 1993, where **LATEX EXAMINATION GLOVES POWDERED** are classified as Class I, declare:

That we are manufacturing for our EU Representative:

and we do hereby declare that the above mentioned product, described under the Technical Documentation, do conform to the requirements set out in Annex VII of Council Directive 93/42/EEC, together with the European Normative 455 /1/2/3 ,and EN 420, according with so that our product are allowed to bear the "CE" mark.

June 12th , 2.020

Yours Sincerely,



P. Selvan
Production / Quality Manager