

EU Type-Examination Certificate

Certificate number: 2777/14815-01/E00-00

This EU Type-Examination Certificate covers the following product group(s) supported by testing to the relevant standards/technical specifications and examination of the technical file documentation:

Following the EU Type-Examination this product group has been shown to satisfy the applicable essential health and safety requirements of Annex II of the PPE Regulation (EU) 2016/425 as a Category III product.

Product reference:

Blue NGPF7000-7005
 Black BNPF7000-7005
 White WNPF7000-7005
 Violet NGPF7000-V1-7005-V1

Description:

Disposable Nitrile gloves(Non-sterile).

Sizes:

6-11(XS-XXL)

Classification:

EN ISO 374-1:2016+A1:2018/Type B

Sodium hydroxide 40%(K)
 Hydrogen peroxide 30%(P)
 Formaldehyde 37%(T)

Level

6
 2
 5

EN ISO 374-4:2019 Degradation %

-68.1
 30.5
 9.5

EN ISO 374-15:2016

Protection against bacteria and fungi
 Protection against viruses

Level

Pass
 Pass

Standards/Technical specifications applied:

EN ISO 21420:2020; EN ISO 374-1:2016+A1:2018; EN ISO 374-5:2016

Technical reports/Approval documents:

SATRA: CHT0296241/2012, CHM0298100/2020/EN/A, CHM0298100/2020/EN/B
 SGS: CH:TX:1142011147, CH:TX:1142011145-1, CH:TX:1142011148
 TUV: 7191234075-CHM20-02-TSL, 7191235025-EEC20-WBH_CR1. 721652920

Signed on behalf of SATRA:

Ting Huang

Ting Huang

Jacque Glasspool

Jacque Glasspool

Date first issued: 20/07/2020

Date of issue: 20/07/2020

Expiry date: 20/07/2025



Packaging Dimensions (mm) & Loading Capabilities (Each product may differ from others)	L x W x H	Method 1 (1000 PCs/CTN)
	Dispenser	230x125x60mm
	Trans. CTN.	315x258x250mm
	Capability for 20 GP	1500 CTNs
	Capability for 40 GP	3100 CTNs
	Capability for 40 HQ	3450 CTNs