



PRICE LIST

We are glad to present you **PPS Gloves Safe Care**. Disposable nitrile gloves designed, tested, and certified in the European Union.

Product specifications

Type: disposable, non-sterile, powder-free nitrile gloves

Material: Nitrile butadiene rubber (NBR)

Colour: light blue

Glove surface: textured fingertip

Weight and thickness: 3.5G / palm - 0.07 mm ; finger - 0.09 mm

AQL: 1,5

Shelf life: 5years (in the original packaging in accordance with ISO 2230)

Package: 100 pos. gloves x 10 boxes x 1 cartón

Size /overall length: XS, S, M, L, XL/ 240 mm

Conformity to standards: EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 1186, EN 1149:2, ASTM D6319, ASTM D7329

Intended Purpose: Medical device class I as set out in MD Regulation (EU) 2017/745, Single-use protective glove - PPE category III (protection against Chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425, Suitable for food contact (as per Regulation (EC) 1935/2004prTT7 Food Contact Materials)



Packaging & Logistic information

	Dispenser box	Transport cartón	Pallet
Packaging	100 pcs. of gloves per box	10 boxes per cartón	96 cartons per pallet 8 layers of 12 cartons
Dimensions	20,5 x 11 x 6 CM	32 x 23 x 21 CM	120 x 80 x 180 CM

Size	Product code	EAN dispenser box	EAN cartón	Cartón gross weight
xs	PPS-SC-01	3 800 501 874127	3 800 501 874 080	3.6 kg.
s	PPS-SC-02	3 800 501 874158	3 800 501 874 073	3.9 kg.
M	PPS-SC-03	3 800 501 874165	3 800 501 874 066	4.2 kg.
L	PPS-SC-04	3 800 501 874 004	3 800 501 874 059	4.5 kg.
XL	PPS-SC-05	3 800 501 874 042	3 800 501 874 035	4.8 kg.
Label	EN - BG - DE - FR - ES			
HS Code	40151900			